

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF nPoint
 ANY REPRODUCTION IN PART OR WHOLE WITHOUT WRITTEN PERMISSION IS PROHIBITED.

4X $\phi 3.2 \nabla$ THRU ALL
 $\square \phi 6.0 \nabla 10.50$
 STAGE MOUNTING HOLES

0.30
 MOUNTING SURFACE

$\square 140.0$ $\square 130.0$ $\square 68.0$


3X $\square 5.0$

$\square 68.0$
 APERTURE $\square 94.0$

3X M2.5X0.45 - 6H ∇ 5.00
 SPECIMEN MOUNTING HOLES

34.0
 CABLE EXIT

18.7
 27.7
 38.3

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS TOLERANCES ARE: .XXX \pm .02 \pm .5° .XX \pm .08 .X \pm .15 DO NOT SCALE DRAWING CAD FILE: 3715135 INSTALLATION DRAWING.slddrw	 nPoint, Inc. 3030 Laura Lane Suite 100 Middleton, WI 53562 Ph: (608) 824-1770 Fx: (608) 824-1774 www.npoint.com	PART DESCRIPTION NPXY100Z100-135 3-AXIS SCANNER XY 100um Z 100um TRAVEL	
		PART NUMBER 3715135	REV. C
SIZE B		SCALE 1:1	SHEET 1 OF 1