



NOTES:

1. SEE FRONT VIEW FOR DEPTH DIMENSIONS
2. SEE FRONT VIEW FOR DEPTH DIMENSIONS
3. SEE FRONT VIEW FOR DEPTH DIMENSIONS
4. SEE FRONT VIEW FOR DEPTH DIMENSIONS
5. ~5N CLAMPING FORCE FOR FRAME # 00-55-500-1000

UNLESS OTHERWISE NOTED, ALL UNITS ARE IN MILLIMETERS

THIRD ANGLE PROJECTION

nPoint, Inc.
3030 Laura Ln. Ste. 100
Middleton, WI 53562

Phone: (608) 824-1770
Fax: (608) 824-1774
www.npoint.com

DESCRIPTION		
NPZ300-222		
1-AXIS SCANNER 300um		
86.5x128.5 APERTURE		
NUMBER	REV.	SHEET 1 OF 1
3715222	A	