

SAMPLE MOUNTING HOLES

8X M2.5X0.45 - 6H  $\nabla$  5

STAGE MOUNTING HOLES

3X M4X7 THRU  
 $\square$   $\phi$  6.5  $\nabla$  19.20

$\phi$  40  
 APERTURE

CH2

CH1

4X 3 X 45°


NOTE  $\triangle$  X

1. THE STAGE MOVES 250 $\mu$ m OUTSIDE THE 75mm FRAME DURING OPERATION IN THE X AND Y AXIS. CARE MUST BE TAKEN DURING MOUNTING OF THE STAGE TO ENSURE ADEQUATE CLEARANCE FOR THIS MOVEMENT.
2. DO NOT CLAMP THE STAGE BODY.

CABLE EXIT

16

UNLESS OTHERWISE NOTED, ALL UNITS ARE IN MILLIMETERS



THIRD ANGLE PROJECTION



nPoint, Inc.  
 3030 Laura Ln. Ste. 100  
 Middleton, WI 53562  
 Phone: (608) 824-1770  
 Fax: (608) 824-1774  
 www.npoint.com

APPROVALS		DATE	
DRAWN: DB		12/04/2014	
DESCRIPTION			
NPXY250 - 405			
NUMBER	REV.	SHEET 1 OF 1	
3715405	A		

